



Human SGSM3 blocking peptide (CDBP2583)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-RUTBC3 antibody
Antigen Description	SGSM3 (small G protein signaling modulator 3) is a protein-coding gene. Diseases associated with SGSM3 include visual agnosia, and agnosia. GO annotations related to this gene include Rab GTPase activator activity and Rab GTPase binding. An important paralog of this gene is TBC1D2.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	SGSM3 small G protein signaling modulator 3 [Homo sapiens (human)]
Official Symbol	SGSM3
Synonyms	SGSM3; small G protein signaling modulator 3; MAP; RUSC3; RUTBC3; RABGAP5; RabGAP-5; rabGAPLP; merlin binding protein; RUN and SH3 containing 3; merlin-associated protein; RUN and TBC1 domain containing 3; RUN and TBC1 domain-containing protein 3; rab-GTPase-activating protein-like protein; small G protein signaling modulator 3 protein;

Entrez Gene ID	27352
mRNA Refseq	NM_015705.4
Protein Refseq	NP_056520.2
UniProt ID	B9A6J5
Chromosome Location	22q13.1-q13.2
Function	Rab GTPase activator activity; Rab GTPase binding; protein binding;