



Mouse SLFN8 blocking peptide (CDBP2614)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-Schlafen 8 (Mouse) antibody
Antigen Description	The Slfn8 gene is conserved in human, chimpanzee, Rhesus monkey, and rat.
Species	Mouse
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	Slfn8 schlafen 8 [Mus musculus (house mouse)]
Official Symbol	SLFN8
Synonyms	SLFN8; schlafen 8;
Entrez Gene ID	276950
mRNA Refseq	NM_001167743.1
Protein Refseq	NP_001161215.1

UniProt ID	E9Q0R4
Chromosome Location	11 C; 11
Function	molecular_function;