



Human SCR3 blocking peptide (CDBP2623)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-SCR3 antibody
Antigen Description	SCR3 (secernin 3) is a protein-coding gene. GO annotations related to this gene include dipeptidase activity. An important paralog of this gene is SCR2.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	SCR3 secernin 3 [Homo sapiens]
Official Symbol	SCR3
Synonyms	SCR3; secernin 3; secernin-3; FLJ23142; SES3; MGC149597;
Entrez Gene ID	79634

mRNA Refseq	NM_001193528
Protein Refseq	NP_001180457
UniProt ID	Q0VDG4
Chromosome Location	2q31
Function	dipeptidase activity;