



Human SESTD1 blocking peptide (CDBP2651)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-SESTD1 antibody
Antigen Description	SESTD1 (SEC14 and spectrin domains 1) is a protein-coding gene. GO annotations related to this gene include phosphatidylinositol-4,5-bisphosphate binding and phosphatidylinositol-3-phosphate binding.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	SESTD1 SEC14 and spectrin domains 1 [Homo sapiens]
Official Symbol	SESTD1
Synonyms	SESTD1; SEC14 and spectrin domains 1; SEC14 domain and spectrin repeat-containing protein 1; DKFZp434O0515; Solo; protein Solo; huntingtin-interacting protein-like protein;

SOLO;

Entrez Gene ID	91404
mRNA Refseq	NM_178123
Protein Refseq	NP_835224
UniProt ID	Q86VW0
Chromosome Location	2q31.3
Function	NOT 1-phosphatidylinositol binding; phosphatidic acid binding; NOT phosphatidylcholine binding; phosphatidylinositol-3,4-bisphosphate binding; phosphatidylinositol-3,5-bisphosphate binding; phosphatidylinositol-3-phosphate binding; phosphatidylinositol-4,