



Human SH3BP1 blocking peptide (CDBP2662)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-SH3BP1 antibody
Antigen Description	SH3BP1 (SH3-domain binding protein 1) is a protein-coding gene. Diseases associated with SH3BP1 include wiskott-aldrich syndrome. GO annotations related to this gene include SH3 domain binding and GTPase activator activity. An important paralog of this gene is ARHGAP17.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	SH3BP1 SH3-domain binding protein 1 [Homo sapiens]
Official Symbol	SH3BP1
Synonyms	SH3BP1; SH3-domain binding protein 1; SH3 domain-binding protein 1; ARHGAP43; 3BP-1; SH3 domain binding protein 1;
Entrez Gene ID	23616
mRNA Refseq	NM_018957

Protein Refseq	NP_061830
UniProt ID	Q9Y3L3
Chromosome Location	22q13.1
Function	GTPase activator activity; SH3 domain binding;