



Human SHISA9 blocking peptide (CDBP2665)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-SHISA9 antibody
Antigen Description	SHISA9 (shisa family member 9) is a protein-coding gene. Diseases associated with SHISA9 include schizophrenia, and neuronitis. An important paralog of this gene is SHISA6.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	SHISA9 shisa homolog 9 (Xenopus laevis) [Homo sapiens]
Official Symbol	SHISA9
Synonyms	SHISA9; shisa homolog 9 (Xenopus laevis); protein shisa-9; cystine-knot AMPAR modulating protein; CKAMP44; FLJ10015; FLJ37653; FLJ58277; FLJ99079; FLJ99234; DKFZp686D24206; DKFZp686E24206;

Entrez Gene ID	729993
mRNA Refseq	NM_001145204
Protein Refseq	NP_001138676
UniProt ID	B4DS77
Chromosome Location	16p13.12