



Human SIDT1 blocking peptide (CDBP2670)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-SIDT1 antibody
Antigen Description	SIDT1 (SID1 transmembrane family, member 1) is a protein-coding gene. GO annotations related to this gene include RNA transmembrane transporter activity. An important paralog of this gene is SIDT2.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	SIDT1 SID1 transmembrane family, member 1 [Homo sapiens]
Official Symbol	SIDT1
Synonyms	SIDT1; SID1 transmembrane family, member 1; SID1 transmembrane family member 1; FLJ20174; SID 1; SID1; SID-1; B830021E24Rik;
Entrez Gene ID	54847
mRNA Refseq	NM_017699

Protein Refseq	NP_060169
UniProt ID	Q9NXL6
Chromosome Location	3q13.31