



# Human SLAIN2 blocking peptide (CDBP2693)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Blocking/Immunizing peptide for anti-SLAIN2 antibody
<b>Antigen Description</b>	SLAIN2 (SLAIN motif family, member 2) is a protein-coding gene. Diseases associated with SLAIN2 include intrahepatic cholangiocarcinoma, and cholangiocarcinoma. GO annotations related to this gene include protein binding. An important paralog of this gene is SLAIN1.
<b>Species</b>	Human
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Apuri, BL, ELISA
<b>Format</b>	Lyophilized powder
<b>Size</b>	100 µg
<b>Preservative</b>	None
<b>Storage</b>	Shipped at ambient temperature, store at -20°C.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">SLAIN2 SLAIN motif family, member 2 [ Homo sapiens ]</a>
<b>Official Symbol</b>	SLAIN2
<b>Synonyms</b>	SLAIN2; SLAIN motif family, member 2; KIAA1458; SLAIN motif-containing protein 2; FLJ21611;
<b>Entrez Gene ID</b>	<a href="#">57606</a>
<b>mRNA Refseq</b>	<a href="#">NM_020846</a>

<b>Protein Refseq</b>	<a href="#">NP_065897</a>
<b>UniProt ID</b>	Q9P270
<b>Chromosome Location</b>	4p12
<b>Function</b>	protein binding;