



Human SNX8 blocking peptide (CDBP2761)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-SNX8 antibody
Antigen Description	SNX8 (sorting nexin 8) is a protein-coding gene. Diseases associated with SNX8 include tetralogy of fallot. GO annotations related to this gene include protein binding and phosphatidylinositol binding.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	SNX8 sorting nexin 8 [Homo sapiens]
Official Symbol	SNX8
Synonyms	SNX8; sorting nexin 8; sorting nexin-8; Mvp1; DKFZp761E1721;
Entrez Gene ID	29886
mRNA Refseq	NM_013321

Protein Refseq	NP_037453
UniProt ID	Q9Y5X2
Chromosome Location	7p22.3
Function	lipid binding; phosphatidylinositol binding; protein binding;