



Human STOX2 blocking peptide (CDBP2855)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-STOX2 antibody
Antigen Description	STOX2 (storkhead box 2) is a protein-coding gene. An important paralog of this gene is STOX1.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	STOX2 storkhead box 2 [Homo sapiens (human)]
Official Symbol	STOX2
Synonyms	STOX2; storkhead box 2; storkhead-box protein 2;
Entrez Gene ID	56977

mRNA Refseq	NM_020225.1
Protein Refseq	NP_064610.1
UniProt ID	Q9P2F5
Chromosome Location	4q35.1
Function	molecular_function;