



Human UBASH3A blocking peptide (CDBP2859)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-STS2/TULA (internal) antibody
Antigen Description	This gene encodes one of two family members belonging to the T-cell ubiquitin ligand (TULA) family. Both family members can negatively regulate T-cell signaling. This family member can facilitate growth factor withdrawal-induced apoptosis in T cells, which may occur via its interaction with AIF, an apoptosis-inducing factor. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Aug 2011]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	UBASH3A ubiquitin associated and SH3 domain containing A [Homo sapiens]
Official Symbol	UBASH3A
Synonyms	UBASH3A; ubiquitin associated and SH3 domain containing A; ubiquitin-associated and SH3

domain-containing protein A; CLIP4; STS 2; TULA; T-cell ubiquitin ligand 1; cbl-interacting protein 4; T-cell ubiquitin ligand protein; suppressor of T-cell receptor signaling 2; STS-2; TULA-1;

Entrez Gene ID	53347
----------------	-----------------------

mRNA Refseq	NM_001001895
-------------	------------------------------

Protein Refseq	NP_001001895
----------------	------------------------------

UniProt ID	P57075
------------	--------

Chromosome Location	21q22.3
---------------------	---------
