



Human SV2A blocking peptide (CDBP2876)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-SV2A antibody
Antigen Description	SV2A (synaptic vesicle glycoprotein 2A) is a protein-coding gene. Diseases associated with SV2A include gastrointestinal stromal tumor, and focal epilepsy, and among its related super-pathways are Non-integrin membrane-ECM interactions. GO annotations related to this gene include transmembrane transporter activity and receptor activity. An important paralog of this gene is SV2B.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	SV2A synaptic vesicle glycoprotein 2A [Homo sapiens]
Official Symbol	SV2A
Synonyms	SV2A; synaptic vesicle glycoprotein 2A; KIAA0736; SV2;
Entrez Gene ID	9900

mRNA Refseq	NM_014849
Protein Refseq	NP_055664
UniProt ID	Q7L0J3
Chromosome Location	1q21.2
Pathway	ECM-receptor interaction, organism-specific biosystem; ECM-receptor interaction, conserved biosystem;
Function	receptor activity; transmembrane transporter activity;