



Human SWAP70 blocking peptide (CDBP2877)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-SWAP70 antibody
Antigen Description	SWAP70 (SWAP switching B-cell complex 70kDa subunit) is a protein-coding gene. Diseases associated with SWAP70 include nodular lymphocyte predominant hodgkin lymphoma, and x-linked hyper igm syndrome. GO annotations related to this gene include phospholipid binding and ATP binding. An important paralog of this gene is DEF6.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	SWAP70 SWAP switching B-cell complex 70kDa subunit [Homo sapiens]
Official Symbol	SWAP70
Synonyms	SWAP70; SWAP switching B-cell complex 70kDa subunit; switch-associated protein 70; KIAA0640; SWAP 70; HSPC321; SWAP-70; FLJ39540;
Entrez Gene ID	23075

mRNA Refseq	NM_015055
Protein Refseq	NP_055870
UniProt ID	Q9UH65
Chromosome Location	11p15
Function	ATP binding; DNA binding; calcium ion binding;