



Human SYNGR4 blocking peptide (CDBP2882)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-Synaptogyrin 4 antibody
Antigen Description	This gene encodes an integral membrane protein. The gene belongs to the synaptogyrin gene family. Like other members of the family the protein contains four transmembrane regions. The exact function of this protein is unclear. [provided by RefSeq, Jul 2008]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	SYNGR4 synaptogyrin 4 [Homo sapiens (human)]
Official Symbol	SYNGR4
Synonyms	SYNGR4; synaptogyrin 4; synaptogyrin-4;
Entrez Gene ID	23546
mRNA Refseq	NM_012451.3

Protein Refseq	NP_036583.2
UniProt ID	O95473
Chromosome Location	19q13.3