



Human SYPL1 blocking peptide (CDBP2894)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-SYPL1 antibody
Antigen Description	SYPL1 (synaptophysin-like 1) is a protein-coding gene. Diseases associated with SYPL1 include alcoholism, and neuronitis. GO annotations related to this gene include transporter activity. An important paralog of this gene is SYNPR.
Species	Human
Conjugate	Unconjugated
Applications	BL, WB
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	SYPL1 synaptophysin-like 1 [Homo sapiens]
Official Symbol	SYPL1
Synonyms	SYPL1; synaptophysin-like 1; synaptophysin like protein , SYPL; synaptophysin-like protein 1; pantophysin; SYPL; H-SP1;

Entrez Gene ID	6856
mRNA Refseq	NM_006754
Protein Refseq	NP_006745
UniProt ID	Q16563
Chromosome Location	7q22.2
Function	transporter activity;