



Human SYPL2 blocking peptide (CDBP2895)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-SYPL2 antibody
Antigen Description	SYPL2 (synaptophysin-like 2) is a protein-coding gene. Diseases associated with SYPL2 include major depressive disorder, and tetanus. GO annotations related to this gene include transporter activity. An important paralog of this gene is SYNPR.
Species	Human
Conjugate	Unconjugated
Applications	BL, WB
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	SYPL2 synaptophysin-like 2 [Homo sapiens]
Official Symbol	SYPL2
Synonyms	SYPL2; synaptophysin-like 2; synaptophysin-like protein 2; Mg29; mitsugumin 29; mitsugumin-29; MG29;

Entrez Gene ID	284612
mRNA Refseq	NM_001040709
Protein Refseq	NP_001035799
UniProt ID	Q5VXT5
Chromosome Location	1p13.3
Function	transporter activity;