



Human TNKS1BP1 blocking peptide (CDBP2898)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-TAB182 antibody
Antigen Description	TNKS1BP1 (tankyrase 1 binding protein 1, 182kDa) is a protein-coding gene. Diseases associated with TNKS1BP1 include prostatitis, and prostate cancer. GO annotations related to this gene include ankyrin binding and enzyme binding.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	TNKS1BP1 tankyrase 1 binding protein 1, 182kDa [Homo sapiens (human)]
Official Symbol	TNKS1BP1
Synonyms	TNKS1BP1; tankyrase 1 binding protein 1, 182kDa; TAB182; 182 kDa tankyrase-1-binding protein; tankyrase 1-binding protein of 182 kDa;
Entrez Gene ID	85456

mRNA Refseq	NM_033396.2
Protein Refseq	NP_203754.2
UniProt ID	Q9C0C2
Chromosome Location	11q12.1
Pathway	Deadenylation of mRNA, organism-specific biosystem; Deadenylation-dependent mRNA decay, organism-specific biosystem; Gene Expression, organism-specific biosystem;
Function	ankyrin binding; enzyme binding;