



Human DCAF12 blocking peptide (CDBP2922)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-TCC52 antibody
Antigen Description	This gene encodes a WD repeat-containing protein that interacts with the COP9 signalosome, a macromolecular complex that interacts with cullin-RING E3 ligases and regulates their activity by hydrolyzing cullin-Nedd8 conjugates. [provided by RefSeq, Jul 2009]
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	DCAF12 DDB1 and CUL4 associated factor 12 [Homo sapiens (human)]
Official Symbol	DCAF12
Synonyms	DCAF12; DDB1 and CUL4 associated factor 12; CT102; TCC52; WDR40A; KIAA1892; DDB1- and CUL4-associated factor 12; WD repeat domain 40A; cancer/testis antigen 102; WD repeat-

containing protein 40A; centrosome-related protein TCC52; testis cancer centrosome-related protein;

Entrez Gene ID	25853
----------------	-----------------------

mRNA Refseq	NM_015397.3
-------------	-----------------------------

Protein Refseq	NP_056212.1
----------------	-----------------------------

UniProt ID	Q5T6F0
------------	--------

Chromosome Location	9p13.3
---------------------	--------
