



Human TEX11 blocking peptide (CDBP2950)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|--|
| Product Overview | Blocking/Immunizing peptide for anti-TEX11 antibody |
| Antigen Description | This gene is X-linked and is expressed in only male germ cells. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008] |
| Species | Human |
| Conjugate | Unconjugated |
| Applications | Apuri, BL, ELISA |
| Format | Lyophilized powder |
| Size | 100 µg |
| Preservative | None |
| Storage | Shipped at ambient temperature, store at -20°C. |

GENE INFORMATION

| | |
|------------------------|--|
| Gene Name | TEX11 testis expressed 11 [Homo sapiens] |
| Official Symbol | TEX11 |
| Synonyms | TEX11; testis expressed 11; testis expressed sequence 11; testis-expressed sequence 11 protein; TGC1; TSGA3; |
| Entrez Gene ID | 56159 |
| mRNA Refseq | NM_001003811 |

| | |
|----------------------------|------------------------------|
| Protein Refseq | NP_001003811 |
| UniProt ID | Q8IYF3 |
| Chromosome Location | Xp11 |
| Function | protein binding; |