



Human TFAP2D blocking peptide (CDBP2953)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-TFAP2D antibody
Antigen Description	TFAP2D (transcription factor AP-2 delta (activating enhancer binding protein 2 delta)) is a protein-coding gene. Diseases associated with TFAP2D include parkinson's disease, and prostatitis. GO annotations related to this gene include DNA binding and sequence-specific DNA binding transcription factor activity. An important paralog of this gene is TFAP2C.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	TFAP2D transcription factor AP-2 delta (activating enhancer binding protein 2 delta) [Homo sapiens]
Official Symbol	TFAP2D
Synonyms	TFAP2D; transcription factor AP-2 delta (activating enhancer binding protein 2 delta); TFAP2BL1, transcription factor AP 2 beta (activating enhancer binding protein 2 beta) like 1; transcription factor AP-2-delta; AP-2 like; AP2-delta; transcription factor AP-2-beta-like 1;

activating enhancer-binding protein 2-delta; activating enhancer binding protein 2 beta-like 1; transcription factor AP-2 beta (activating enhancer binding protein 2 beta)-like 1; transcription factor AP-2 beta (activating enhancer-binding protein 2 beta)-like 1; TFAP2BL1;

Entrez Gene ID	83741
----------------	-----------------------

mRNA Refseq	NM_172238
-------------	---------------------------

Protein Refseq	NP_758438
----------------	---------------------------

UniProt ID	Q7Z6R9
------------	--------

Chromosome Location	6p12.3
---------------------	--------

Function	DNA binding; sequence-specific DNA binding transcription factor activity;
----------	---
