



Human HTATIP2 blocking peptide (CDBP2986)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-TIP30/CC3 antibody
Antigen Description	HTATIP2 (HIV-1 Tat interactive protein 2, 30kDa) is a protein-coding gene. Diseases associated with HTATIP2 include gnathodiaphyseal dysplasia, and hiv-1. GO annotations related to this gene include transcription coactivator activity and oxidoreductase activity.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	HTATIP2 HIV-1 Tat interactive protein 2, 30kDa [Homo sapiens (human)]
Official Symbol	HTATIP2
Synonyms	HTATIP2; HIV-1 Tat interactive protein 2, 30kDa; CC3; TIP30; SDR44U1; oxidoreductase HTATIP2; Tat-interacting protein (30kD); HIV-1 TAT-interactive protein 2; 30 kDa HIV-1 TAT-interacting protein; short chain dehydrogenase/reductase family 44U, member 1;
Entrez Gene ID	10553

mRNA Refseq	NM_001098520.1
Protein Refseq	NP_001091990.1
UniProt ID	Q9BUP3
Chromosome Location	11p15.1
Function	oxidoreductase activity; protein binding; transcription coactivator activity;