



Human PLIN3 blocking peptide (CDBP2987)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-TIP47 antibody
Antigen Description	HTATIP2 (HIV-1 Tat interactive protein 2, 30kDa) is a protein-coding gene. Diseases associated with HTATIP2 include gnathodiaphyseal dysplasia, and hiv-1. GO annotations related to this gene include transcription coactivator activity and oxidoreductase activity.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	PLIN3 perilipin 3 [Homo sapiens (human)]
Official Symbol	PLIN3
Synonyms	PLIN3; perilipin 3; PP17; TIP47; M6PRBP1; perilipin-3; placental protein 17; 47 kDa MPR-binding protein; cargo selection protein TIP47; tail-interacting protein, 47 kD; mannose-6-

phosphate receptor-binding protein 1;

Entrez Gene ID	10226
mRNA Refseq	NM_001164189.1
Protein Refseq	NP_001157661.1
UniProt ID	O60664
Chromosome Location	19p13.3
Function	protein binding;