



Human TIPRL blocking peptide (CDBP2989)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-TIPRL antibody
Antigen Description	TIPRL is an inhibitory regulator of protein phosphatase-2A (PP2A) (see PPP2CA; MIM 176915), PP4 (see PPP4C; MIM 602035), and PP6 (see PPP6C; MIM 612725) (McConnell et al., 2007 [PubMed 17384681]).[supplied by OMIM, Nov 2010]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	TIPRL TIP41, TOR signaling pathway regulator-like (S. cerevisiae) [Homo sapiens]
Official Symbol	TIPRL
Synonyms	TIPRL; TIP41, TOR signaling pathway regulator-like (S. cerevisiae); TIP41, TOR signalling pathway regulator like (S. cerevisiae); TIP41-like protein; dJ69E11.3; MGC3794; TIP41; type 2A-interacting protein; putative MAPK activating protein; putative MAPK-activating protein PM10; TIP; RP1-69E11.1;
Entrez Gene ID	261726

mRNA Refseq	NM_001031800
Protein Refseq	NP_001026970
UniProt ID	O75663
Chromosome Location	1q23.2
Function	protein binding;