



Human TMEM184B blocking peptide (CDBP3000)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-TMEM184B antibody
Antigen Description	TMEM184B (transmembrane protein 184B) is a protein-coding gene. An important paralog of this gene is TMEM184C.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	TMEM184B transmembrane protein 184B [Homo sapiens]
Official Symbol	TMEM184B
Synonyms	TMEM184B; transmembrane protein 184B; C22orf5, chromosome 22 open reading frame 5;

DKFZP586A1024; FM08; HS5O6A; putative MAPK-activating protein FM08; C22orf5;
DKFZp586A1024; DKFZp686N1150; DKFZp686F07174;

Entrez Gene ID	25829
mRNA Refseq	NM_001195071
Protein Refseq	NP_001182000
UniProt ID	Q9Y519
Chromosome Location	22q12