



Human ZRANB1 blocking peptide (CDBP3031)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-TRABID/ZRANB1 antibody
Antigen Description	ZRANB1 (zinc finger, RAN-binding domain containing 1) is a protein-coding gene. GO annotations related to this gene include K63-linked polyubiquitin binding and ubiquitin thioesterase activity. An important paralog of this gene is TNFAIP3.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	ZRANB1 zinc finger, RAN-binding domain containing 1 [Homo sapiens]
Official Symbol	ZRANB1
Synonyms	ZRANB1; zinc finger, RAN-binding domain containing 1; ubiquitin thioesterase ZRANB1; TRABID; hTrabid; TRAF-binding protein domain; TRAF-binding domain-containing protein; zinc finger Ran-binding domain-containing protein 1; DKFZp762P2216;
Entrez Gene ID	54764

mRNA Refseq	NM_017580
Protein Refseq	NP_060050
UniProt ID	Q9UGI0
Chromosome Location	10q26.12
Function	K63-linked polyubiquitin binding; cysteine-type peptidase activity; metal ion binding; peptidase activity; protein binding; ubiquitin thiolesterase activity; ubiquitin-specific protease activity; zinc ion binding;