



Human TREX2 blocking peptide (CDBP3047)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-TREX2 antibody
Antigen Description	This gene encodes a nuclear protein with 3' to 5' exonuclease activity. The encoded protein participates in double-stranded DNA break repair, and may interact with DNA polymerase delta. [provided by RefSeq, Nov 2012]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	TREX2 three prime repair exonuclease 2 [Homo sapiens]
Official Symbol	TREX2
Synonyms	TREX2; three prime repair exonuclease 2; 3-5 exonuclease TREX2 long form;
Entrez Gene ID	11219
mRNA Refseq	NM_080701

Protein Refseq	NP_542432
UniProt ID	Q9BQ50
Chromosome Location	Xq28
Function	3-5-exodeoxyribonuclease activity; exodeoxyribonuclease III activity; exonuclease activity; hydrolase activity; nucleic acid binding;