



Human TRIM71 blocking peptide (CDBP3067)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-TRIM71 (aa409-419) antibody
Antigen Description	TRIM71 (tripartite motif containing 71, E3 ubiquitin protein ligase) is a protein-coding gene. Diseases associated with TRIM71 include otitis media, and cholesteatoma. GO annotations related to this gene include ubiquitin-protein ligase activity and miRNA binding. An important paralog of this gene is TRIM3.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	TRIM71 tripartite motif containing 71, E3 ubiquitin protein ligase [Homo sapiens]
Official Symbol	TRIM71
Synonyms	TRIM71; tripartite motif containing 71, E3 ubiquitin protein ligase; tripartite motif containing 71 , tripartite motif containing 71; tripartite motif-containing protein 71; LIN 41; LIN41; protein lin-41 homolog; abnormal cell LiNeage LIN-41; homolog of C. elegans Lin-41; tripartite motif-containing 71; LIN-41;

Entrez Gene ID	131405
mRNA Refseq	NM_001039111
Protein Refseq	NP_001034200
UniProt ID	Q2Q1W2
Chromosome Location	3p22.3
Function	metal ion binding; zinc ion binding;