



Human TRIM71 blocking peptide (CDBP3068)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-TRIM71 antibody
Antigen Description	TRIM71 (tripartite motif containing 71, E3 ubiquitin protein ligase) is a protein-coding gene. Diseases associated with TRIM71 include otitis media, and cholesteatoma. GO annotations related to this gene include ubiquitin-protein ligase activity and miRNA binding. An important paralog of this gene is TRIM3.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	TRIM71 tripartite motif containing 71, E3 ubiquitin protein ligase [Homo sapiens]
Official Symbol	TRIM71
Synonyms	TRIM71; tripartite motif containing 71, E3 ubiquitin protein ligase; tripartite motif containing 71 ,

tripartite motif containing 71; tripartite motif-containing protein 71; LIN 41; LIN41; protein lin-41 homolog; abnormal cell LIneage LIN-41; homolog of C. elegans Lin-41; tripartite motif-containing 71; LIN-41;

Entrez Gene ID	131405
----------------	------------------------

mRNA Refseq	NM_001039111
-------------	------------------------------

Protein Refseq	NP_001034200
----------------	------------------------------

UniProt ID	Q2Q1W2
------------	--------

Chromosome Location	3p22.3
---------------------	--------

Function	metal ion binding; zinc ion binding;
----------	--------------------------------------
