



Human TTC5 blocking peptide (CDBP3097)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-TTC5 antibody
Antigen Description	TTC5 (tetratricopeptide repeat domain 5) is a protein-coding gene. GO annotations related to this gene include chromatin binding and DNA binding.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	TTC5 tetratricopeptide repeat domain 5 [Homo sapiens]
Official Symbol	TTC5
Synonyms	TTC5; tetratricopeptide repeat domain 5; tetratricopeptide repeat protein 5; Strap; TPR repeat protein 5;

Entrez Gene ID	91875
mRNA Refseq	NM_138376
Protein Refseq	NP_612385
UniProt ID	Q8N0Z6
Chromosome Location	14q11.2
Pathway	p53 pathway, organism-specific biosystem;
Function	binding;