



Human UCHL1 blocking peptide (CDBP3120)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-UCHL1 (aa 211-221) antibody
Antigen Description	The protein encoded by this gene belongs to the peptidase C12 family. This enzyme is a thiol protease that hydrolyzes a peptide bond at the C-terminal glycine of ubiquitin. This gene is specifically expressed in the neurons and in cells of the diffuse neuroendocrine system. Mutations in this gene may be associated with Parkinson disease.(provided by RefSeq, Sep 2009)
Species	Human
Conjugate	Unconjugated
Applications	Ahuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	uchl1 ubiquitin carboxyl-terminal esterase L1 (ubiquitin thiolesterase) [Xenopus (Silurana) tropicalis]
Official Symbol	UCHL1
Synonyms	UCHL1; ubiquitin carboxyl-terminal esterase L1 (ubiquitin thiolesterase); ubiquitin carboxyl-terminal esterase L1; park5; pgp9.5; uch-1;

Entrez Gene ID [496780](#)

mRNA Refseq [NM_001011321](#)

Protein Refseq [NP_001011321](#)
