



Human UNC45B blocking peptide (CDBP3129)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-UNC45B antibody
Antigen Description	UNC45B plays a role in myoblast fusion and sarcomere organization (Price et al., 2002 [PubMed 12356907]).[supplied by OMIM, Mar 2008]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	UNC45B unc-45 homolog B (C. elegans) [Homo sapiens]
Official Symbol	UNC45B
Synonyms	UNC45B; unc-45 homolog B (C. elegans); cardiomyopathy associated 4 , CMYA4; protein unc-45 homolog B; UNC45; UNC-45B; striated muscle UNC45; cardiomyopathy associated 4; CMYA4; SMUNC45; FLJ38610; MGC119540; MGC119541;
Entrez Gene ID	146862
mRNA Refseq	NM_001033576

Protein Refseq	NP_001028748
UniProt ID	Q8IWX7
Chromosome Location	17q12
Function	binding;