



Human USP19 blocking peptide (CDBP3145)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-USP19 antibody
Antigen Description	USP19 (ubiquitin specific peptidase 19) is a protein-coding gene. Diseases associated with USP19 include hypoxia. GO annotations related to this gene include ubiquitin-specific protease activity and ubiquitin thiolesterase activity. An important paralog of this gene is USP45.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	USP19 ubiquitin specific peptidase 19 [Homo sapiens]
Official Symbol	USP19
Synonyms	USP19; ubiquitin specific peptidase 19; ubiquitin specific protease 19; ubiquitin carboxyl-terminal hydrolase 19; KIAA0891; ZMYND9; ubiquitin thioesterase 19; deubiquitinating enzyme 19; ubiquitin thiolesterase 19; ubiquitin-specific-processing protease 19; zinc finger MYND domain-containing protein 9;
Entrez Gene ID	10869

mRNA Refseq	NM_001199160
Protein Refseq	NP_001186089
UniProt ID	O94966
Chromosome Location	3p21.31
Function	cysteine-type peptidase activity; metal ion binding; peptidase activity; protein binding; ubiquitin thiolesterase activity; ubiquitin-specific protease activity; zinc ion binding;