



Human USP3 blocking peptide (CDBP3151)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-USP3 antibody
Antigen Description	USP3 (ubiquitin specific peptidase 3) is a protein-coding gene. GO annotations related to this gene include ubiquitin-specific protease activity and ubiquitin thiolesterase activity. An important paralog of this gene is USP44.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	USP3 ubiquitin specific peptidase 3 [Homo sapiens]
Official Symbol	USP3
Synonyms	USP3; ubiquitin specific peptidase 3; ubiquitin specific protease 3; ubiquitin carboxyl-terminal hydrolase 3; ubiquitin thioesterase 3; deubiquitinating enzyme 3; ubiquitin thiolesterase 3; ubiquitin-specific-processing protease 3; UBP; SIH003; MGC129878; MGC129879;
Entrez Gene ID	9960

mRNA Refseq	NM_001256702
Protein Refseq	NP_001243631
UniProt ID	Q9Y6I4
Chromosome Location	15q22.3
Function	cysteine-type peptidase activity; histone binding; metal ion binding; peptidase activity; ubiquitin thiolesterase activity; ubiquitin thiolesterase activity; ubiquitin-specific protease activity; zinc ion binding;