



Human USP43 blocking peptide (CDBP3154)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-USP43 antibody
Antigen Description	USP43 (ubiquitin specific peptidase 43) is a protein-coding gene. GO annotations related to this gene include cysteine-type peptidase activity and ubiquitin thiolesterase activity. An important paralog of this gene is USP45.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	USP43 ubiquitin specific peptidase 43 [Homo sapiens]
Official Symbol	USP43
Synonyms	USP43; ubiquitin specific peptidase 43; ubiquitin specific protease 43; ubiquitin carboxyl-terminal hydrolase 43; FLJ30626; ubiquitin thioesterase 43; deubiquitinating enzyme 43; ubiquitin thiolesterase 43; ubiquitin-specific-processing protease 43;
Entrez Gene ID	124739

mRNA Refseq	NM_153210
Protein Refseq	NP_694942
UniProt ID	Q70EL4
Chromosome Location	17p12
Function	cysteine-type peptidase activity; peptidase activity; ubiquitin thiolesterase activity;