



Human GC blocking peptide (CDBP3161)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-VCBP/GC (aa 18 - 30) antibody
Antigen Description	The protein encoded by this gene belongs to the albumin gene family. It is a multifunctional protein found in plasma, ascitic fluid, cerebrospinal fluid and on the surface of many cell types. It binds to vitamin D and its plasma metabolites and transports them to target tissues. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	GC group-specific component (vitamin D binding protein) [Homo sapiens]
Official Symbol	GC
Synonyms	GC; group-specific component (vitamin D binding protein); vitamin D-binding protein; DBP; hDBP; VDBP; VDB; gc-globulin; vitamin D-binding alpha-globulin; GRD3; VDBG; DBP/GC;
Entrez Gene ID	2638

mRNA Refseq	NM_000583
Protein Refseq	NP_000574
UniProt ID	P02774
Chromosome Location	4q12-q13
Pathway	Metabolism, organism-specific biosystem; Metabolism of lipids and lipoproteins, organism-specific biosystem; Metabolism of steroid hormones and vitamins A and D, organism-specific biosystem; Vitamin D (calciferol) metabolism, organism-specific biosystem;
Function	actin binding; vitamin D binding; vitamin transporter activity;