



Human WBP2NL blocking peptide (CDBP3189)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-WBP2NL/PAWP antibody
Antigen Description	WBP2NL is a sperm-specific WW domain-binding protein that promotes meiotic resumption and pronuclear development during oocyte fertilization (Wu et al., 2007 [PubMed 17289678]).[supplied by OMIM, Mar 2008]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	WBP2NL WBP2 N-terminal like [Homo sapiens]
Official Symbol	WBP2NL
Synonyms	WBP2NL; WBP2 N-terminal like; postacrosomal sheath WW domain-binding protein; FLJ26145; MGC26816; PAWP; postacrosomal sheath WW domain binding protein; WW domain-binding protein 2-like;
Entrez Gene ID	164684

mRNA Refseq	NM_152613
Protein Refseq	NP_689826
UniProt ID	Q6ICG8
Chromosome Location	22q13.2
Function	WW domain binding;