



Human ZCCHC11 blocking peptide (CDBP3231)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-ZCCHC11/PAPD3 antibody
Antigen Description	ZCCHC11 is an RNA uridylyltransferase (EC 2.7.7.52) that uses UTP to add uridines to the 3-prime end of substrate RNA molecules (Jones et al., 2009 (PubMed 19701194)).(supplied by OMIM, Jan 2011)
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	ZCCHC11 zinc finger, CCHC domain containing 11 [Homo sapiens]
Official Symbol	ZCCHC11
Synonyms	ZCCHC11; zinc finger, CCHC domain containing 11; terminal uridylyltransferase 4; KIAA0191; PAPD3; TUTase4; TUTase 4; PAP associated domain containing 3; zinc finger CCHC domain-containing protein 11; TUT4; FLJ42878; DKFZp779C1943;
Entrez Gene ID	23318

mRNA Refseq	NM_001009881
Protein Refseq	NP_001009881
UniProt ID	Q5TAX3
Chromosome Location	1p32.3
Function	RNA uridylyltransferase activity; metal ion binding; nucleic acid binding; nucleotidyltransferase activity; protein binding; transferase activity; zinc ion binding;