



Human ZNF217 blocking peptide (CDBP3246)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-ZNF217 antibody
Antigen Description	ZNF217 (zinc finger protein 217) is a protein-coding gene. Diseases associated with ZNF217 include cystadenocarcinoma, and serous cystadenocarcinoma. GO annotations related to this gene include transcription regulatory region DNA binding and sequence-specific DNA binding transcription factor activity. An important paralog of this gene is ZNF219.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	ZNF217 zinc finger protein 217 [Homo sapiens]
Official Symbol	ZNF217
Synonyms	ZNF217; zinc finger protein 217; ZABC1;
Entrez Gene ID	7764
mRNA Refseq	NM_006526

Protein Refseq	NP_006517
UniProt ID	O75362
Chromosome Location	20q13.2
Function	DNA binding; metal ion binding; protein binding; sequence-specific DNA binding transcription factor activity; transcription regulatory region DNA binding; zinc ion binding;