



Human ZFP112 blocking peptide (CDBP3247)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-ZNF228/ZFP112 antibody
Antigen Description	ZNF112 (zinc finger protein 112) is a protein-coding gene. Among its related super-pathways are Gene Expression.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	ZFP112 zinc finger protein 112 homolog (mouse) [Homo sapiens]
Official Symbol	ZFP112
Synonyms	ZFP112; zinc finger protein 112 homolog (mouse); zinc finger protein 228 , ZNF228; zinc finger protein 112 homolog; zfp-112; zinc finger protein 228; ZNF228;
Entrez Gene ID	7771
mRNA Refseq	NM_001083335

Protein Refseq	NP_001076804
UniProt ID	Q9UJU3
Chromosome Location	19q13.2
Pathway	Gene Expression, organism-specific biosystem; Generic Transcription Pathway, organism-specific biosystem;
Function	DNA binding; metal ion binding; zinc ion binding;