



AFAP1L2 blocking peptide (CDBP5024)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	AFAP1L2 actin filament associated protein 1-like 2 [Homo sapiens (human)]
Official Symbol	AFAP1L2
Synonyms	AFAP1L2; actin filament associated protein 1-like 2; XB130; KIAA1914; CTB-1144G6.4; actin filament-associated protein 1-like 2; AFAP1-like protein 2
Entrez Gene ID	84632
mRNA Refseq	NM_001001936
Protein Refseq	NP_001001936
UniProt ID	Q8N4X5
Function	SH2 domain binding; SH3 domain binding; protein tyrosine kinase activator activity