



AKIRIN2 blocking peptide (CDBP5043)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	AKIRIN2 akirin 2 [Homo sapiens (human)]
Official Symbol	AKIRIN2
Synonyms	AKIRIN2; akirin 2; FBI1; C6orf166; dJ486L4.2; akirin-2; fourteen-three-three beta interactant 1
Entrez Gene ID	55122
mRNA Refseq	NM_018064
Protein Refseq	NP_060534
UniProt ID	Q53H80
Function	enzyme binding