



ANKLE1 blocking peptide (CDBP5059)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ANKLE1 ankyrin repeat and LEM domain containing 1 [Homo sapiens (human)]
Official Symbol	ANKLE1
Synonyms	ANKLE1; ankyrin repeat and LEM domain containing 1; LEM3; LEMD6; ANKRD41; ankyrin repeat and LEM domain-containing protein 1; LEM domain containing 6; LEM-domain containing 3; ankyrin repeat domain-containing protein 41
Entrez Gene ID	126549
mRNA Refseq	NM_001278443
Protein Refseq	NP_001265372
UniProt ID	Q8NAG6