



APOA4 blocking peptide (CDBP5089)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	APOA4 apolipoprotein A-IV [Gallus gallus (chicken)]
Official Symbol	APOA4
Synonyms	APOA4; apolipoprotein A-IV; APOAIV
Entrez Gene ID	395780
mRNA Refseq	NM_204938
Protein Refseq	NP_990269