



APOA5 blocking peptide (CDBP5090)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	APOA5 apolipoprotein A-V [Gallus gallus (chicken)]
Official Symbol	APOA5
Synonyms	APOA5; apolipoprotein A-V
Entrez Gene ID	419802
Pathway	PPAR signaling pathway