



ARHGAP18 blocking peptide (CDBP5100)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ARHGAP18 belongs to a family of Rho (see MIM 165390) GTPase-activating proteins that modulate cell signaling (Potkin et al., 2009 [PubMed 19065146]).[supplied by OMIM, Apr 2010]
Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ARHGAP18 Rho GTPase activating protein 18 [Homo sapiens (human)]
Official Symbol	ARHGAP18
Synonyms	ARHGAP18; Rho GTPase activating protein 18; SENEX; MacGAP; bA307O14.2; rho GTPase-activating protein 18; rho-type GTPase-activating protein 18
Entrez Gene ID	93663
mRNA Refseq	NM_033515
Protein Refseq	NP_277050

UniProt ID	Q8N392
Pathway	Rho GTPase cycle; Signal Transduction; Signaling by Rho GTPases
Function	GTPase activator activity