



ARHGAP39 blocking peptide (CDBP5101)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ARHGAP39 Rho GTPase activating protein 39 [Homo sapiens (human)]
Official Symbol	ARHGAP39
Synonyms	ARHGAP39; Rho GTPase activating protein 39; CrGAP; Vilse; rho GTPase-activating protein 39; crossGAP homolog; RhoGAP93B homolog
Entrez Gene ID	80728
mRNA Refseq	NM_025251
Protein Refseq	NP_079527
UniProt ID	Q9C0H5
Pathway	Axon guidance; Developmental Biology; Inactivation of Cdc42 and Rac; Rho GTPase cycle;

