



ARL15 blocking peptide (CDBP5102)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ARL15 ADP-ribosylation factor-like 15 [Homo sapiens (human)]
Official Symbol	ARL15
Synonyms	ARL15; ADP-ribosylation factor-like 15; ARFRP2; ADP-ribosylation factor-like protein 15; ARF-related protein 2; ADP-ribosylation factor related protein 2; ADP-ribosylation factor-related protein 2
Entrez Gene ID	54622
mRNA Refseq	NM_019087
Protein Refseq	NP_061960
UniProt ID	Q9NXU5

FunctionGTP binding; molecular_function
