



ARL6IP1 blocking peptide (CDBP5103)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ARL6IP1 ADP-ribosylation factor-like 6 interacting protein 1 [Homo sapiens (human)]
Official Symbol	ARL6IP1
Synonyms	ARL6IP1; ADP-ribosylation factor-like 6 interacting protein 1; AIP1; ARMER; SPG61; ARL6IP; ADP-ribosylation factor-like protein 6-interacting protein 1; aip-1; ARL-6-interacting protein 1; apoptotic regulator in the membrane of the endoplasmic reticulum
Entrez Gene ID	23204
mRNA Refseq	NM_015161
Protein Refseq	NP_055976
UniProt ID	Q15041

Functionprotein binding
