



MAP3K5 blocking peptide (CDBP5111)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Mitogen-activated protein kinase (MAPK) signaling cascades include MAPK or extracellular signal-regulated kinase (ERK), MAPK kinase (MKK or MEK), and MAPK kinase kinase (MAPKKK or MEKK). MAPKK kinase/MEKK phosphorylates and activates its downstream protein kinase, MAPK kinase/MEK, which in turn activates MAPK. The kinases of these signaling cascades are highly conserved, and homologs exist in yeast, Drosophila, and mammalian cells. MAPKKK5 contains 1,374 amino acids with all 11 kinase subdomains. Northern blot analysis shows that MAPKKK5 transcript is abundantly expressed in human heart and pancreas. The MAPKKK5 protein phosphorylates and activates MKK4 (aliases SERK1, MAPKK4) in vitro, and activates c-Jun N-terminal kinase (JNK)/stress-activated protein kinase (SAPK) during transient expression in COS and 293 cells; MAPKKK5 does not activate MAPK/ERK. [provided by RefSeq, Jul 2008]
----------------------------	---

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	MAP3K5 mitogen-activated protein kinase kinase kinase 5 [Homo sapiens (human)]
Official Symbol	MAP3K5

Synonyms	MAP3K5; mitogen-activated protein kinase kinase kinase 5; ASK1; MEKK5; MAPKKK5; ASK-1; MEKK 5; MEK kinase 5; MAP/ERK kinase kinase 5; MAPK/ERK kinase kinase 5; apoptosis signal regulating kinase 1; apoptosis signal-regulating kinase 1
Entrez Gene ID	4217
mRNA Refseq	NM_005923
Protein Refseq	NP_005914
UniProt ID	Q99683
Pathway	Amyotrophic lateral sclerosis (ALS); Apoptosis Modulation and Signaling; Cellular Senescence; Cellular responses to stress; Class I PI3K signaling events mediated by Akt; Corticotropin-releasing hormone; HIV-1 Nef: Negative effector of Fas and TNF-alpha; Insulin Signaling
Function	ATP binding; MAP kinase kinase kinase activity; cysteine-type endopeptidase activator activity involved in apoptotic process; magnesium ion binding; protein binding; protein homodimerization activity; protein kinase activity; protein kinase binding; protein phosphatase binding