



ATG101 blocking peptide (CDBP5117)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Used as a blocking peptide in immunoblotting applications.
Format	Liquid
Concentration	200 µg/mL
Size	0.05 mg
Preservative	None
Storage	-20°C

GENE INFORMATION

Gene Name	ATG101 autophagy related 101 [Homo sapiens (human)]
Official Symbol	ATG101
Synonyms	ATG101; autophagy related 101; C12orf44; autophagy-related protein 101; Atg13-interacting protein
Entrez Gene ID	60673
mRNA Refseq	NM_001098673
Protein Refseq	NP_001092143
UniProt ID	Q9BSB4
Function	identical protein binding; protein binding; protein complex binding